# LIFEPAK® CR2

# defibrillator cellular

## **Features**

- · Easy to use
- Cellular connectivity
- Self-monitoring
- Location services
- cprINSIGHT® analysis technology
- Durability

# You can't predict when an emergency strikes, but you can be prepared.

Sudden cardiac arrest can happen to anyone, anywhere. A victim's chance of survival dramatically decreases for every minute without treatment<sup>1</sup>. A nationwide study of law enforcement agencies proved first responder automated external defibrillator (AED) programs can reduce the time to defibrillation and increase a victim's chance of survival from out-of-hospital sudden cardiac arrests<sup>2</sup>.

# Metronome and CPR coaching

Sets an effective pace and audibly guides users.

#### • ClearVoice™ technology

Detects background noise and adjust tones and voice prompts to ensure they can be heard clearly in noisy environments.

#### • QUIK-STEP® electrodes

Peel directly off the base for faster placement.

## • LIFEPAK TOUGH<sup>TM</sup>

IP55 rating for challenging environments.

### Child mode

Child mode delivers lower energy levels appropriate for young children without having to change electrodes.

#### Cellular connectivity

Device will use a 4G cellular connection and/or Wi-Fi to communicate with LIFELINK central  $^{\text{TM}}$ .

#### • LIFELINKcentral AED Program manager

Monitor AED programs by tracking AED status, sending patient data to hospitals, and giving device specific location information.









#### References:

1 Graham R, McCoy M, Schultz A. Strategies to Improve Cardiac Arrest Survival, A Time to Act. Institute of Medicine Report, 2015. 2. Husain S, Eisenberg M, Police AED programs: a systematic review and meta-analysis, Resuscitation. 2013; 84: 1184–1191.

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate Stryker's product offerings. A healthcare professional must always refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: ClearVoice, cprINSIGHT, LIFELINKcentral, LIFEPAK, LIFEKPAK Tough, QUIK-STEP, Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.



Manufactured by: Physio-Control, Inc. 11811 Willows Road NE Redmond, WA, 98052 U.S.A. Toll free 800 442 1142 stryker.com/emergencycare